



THE FUSION[®] OA BRACE

FUSION OA Frame Design

Sleek, contoured design. Exceptionally lightweight material delivers enhanced strength and allows unhindered mobility.

Cushioned Straps

Easily adjustable/removable straps ensure a proper fit, and the hook and loop feature provides dependable fixation. Additional strap padding is durable, comfortable and can be removed for easy cleaning.

Convenient Thumbwheel Dial*

Exclusive thumbwheel dial hinge makes fine-tuning the level of off-loading easy.

ProForm™ Technology

Medial side structure utilizes Zytel™, a strong yet flexible molded polymer that conforms to the leg's shape to provide a contoured fit and superior protection.

PivotPoint Strap Tabs

Pivoting strap tabs are constructed with Zytel for added strength and durability, which also allows the straps to contour to the leg for superior comfort during motion.

Soft, Contoured Frame Pads

Contoured frame padding conforms to the leg, providing increased comfort.





The exclusive thumbwheel dial hinge makes fine-tuning FUSION OA'S level of off-loading easy. It requires no tools, and offers patient-friendly adjustability in .5mm increments.

Osteoarthritis, as a degenerative disease, is a problem affecting both females and males. Wearing an unloading knee brace may improve gait symmetry, decrease symptomatic pain, and increase activity for patients who have unicompartmental osteoarthritis of the knee or other varus knee deformities.

Osteoarthritis bracing can prolong the interval between the onset of symptoms and prospective surgery, and may slow the degenerative process.

FUSION OA Sizing Chart

Size #	Size	6" Below Mid Patella	6" Above Mid Patella
X=1	XS	10 - 12" (25 - 30 cm)	13 - 15" (33 - 39 cm)
X=2	S	12 - 14" (30 - 36 cm)	15.5 - 18" (39 - 46 cm)
X=3	M	14 - 15" (36 - 38 cm)	18 - 19.5" (46 - 50 cm)
X=4	M+	15 - 16" (38 - 41 cm)	19.5 - 21" (50 - 53 cm)
X=5	L	16 - 18" (41 - 46 cm)	21 - 24" (53 - 61 cm)
X=6	XL	18 - 20" (46 - 51 cm)	24 - 27" (61 - 69 cm)
X=7	XXL	20 - 23" (51 - 58 cm)	27 - 30" (69 - 76 cm)

Circumference should be taken 6" (15 cm) above mid patella and 6" (15 cm) below mid patella.



Distributed by:
 Pacific Medical Inc.
 1700 N. Chrisman Road
 Tracy, CA 95304
 Phone: 866.323.8395 Fax: 866.323.8612
 bracing@pacmedical.com



Freedom to perform FUSION OA was created to meet the unique needs of those suffering from osteoarthritis. Designed for maximum pain relief and all-day comfort, FUSION OA also provides the necessary support for ACL, PCL and collateral ligament protection. The revolutionary ProForm technology ensures an intimate, contoured fit, while FUSION OA's convenient thumbwheel dial hinge makes precision off-loading adjustments easy without using a tool. The brace's streamlined, low profile frame also fits discreetly under clothing. Best of all, FUSION OA not only relieves pain, but could also slow osteoarthritic deterioration. Put simply, FUSION OA lets osteoarthritis patients perform at their best.

FUSION OA Braces

Custom Part No.	Description
502900	FUSION OA CUSTOM
502901	FUSION OA CUSTOM, COLOR OPTION
OTS Part No.	
50283X	FUSION OA, OTS, LT, XS - XXL - medial
50284X	FUSION OA, OTS, RT, XS - XXL - medial

Top Shelf Orthopedics
 Tracy, Ca 95304
 Phone: 866.592.0488
 Fax: 209.834.8832

