

The Gold Standard...

Proven design with enhanced features.

- Active brace stability with dynamic strap design provides neutralizing force to the tibia
- Integrated strapping system enhances control and reduces profile
- Lightweight brace accommodates any time wear
- Graduated thigh reduces medial profile and accommodates bilateral brace wear
- Medial tibial contour controls medial and lateral tibial translation

Anterior thigh cuff anchors the femur



Posterior thigh strap pushes the femur forward

Anterior tibial strap provides a neutralizing force on the tibia preventing anterior tibial translation

Posterior cuff anchors the tibia

- Suede moisture wicking liner draws moisture away from skin
- Aerospace grade aluminum provides rigidity for collision sports
- Three configurations and two lengths are available in ACL, Combined Instability, and PCL version with short and standard calf
- Latex free silicone tibial strap pad adds suspension control
- Carry bag included

Control
Fit
Comfort

Active Brace Stability